

# Just Like Money

## Prepaid Calibration Vouchers

**Rainin Prepaid Calibration Vouchers let you purchase pipette service when convenient and redeem them as needed.**

- Good for any brand of pipettes
- Valid for mail-in calibration service only
- Price guaranteed for 2 years
- Free shipping, next business day
  - 4 or more pipettes, request our free return kit at [mt.com/rainin-kit](http://mt.com/rainin-kit)
  - 3 or less, pack them yourself and print a free shipping label at [mt.com/RaininServiceLabels](http://mt.com/RaininServiceLabels)

### Calibration Plans

Standard Plan	Advanced Plan	Accredited Plan
Labs with routine service needs	When prior performance matters	For complete ISO 17025 documentation
Perfect for <b>Academic Research R&amp;D</b>	<b>Pharma Biotech Clinical</b> research labs	<b>GLP/GMP Regulated Labs</b>

See reverse page for detail information



# RAININ

Feel the difference!

## Prepaid Calibration Voucher Details and Ordering Information

STANDARD PLAN – VOUCHER	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
For labs with routine service needs	✓	n/a	4 test weighings @ 10% and 100%	Standard Certificate

- Full preventive maintenance (PM) with seal replacement and leak test
- Adjustment to manufacturer's tolerances
- Calibration certificate with detailed results and historical trending

ORDER No.	DESCRIPTION	PRICE
<b>17014448</b>	Single channel	\$ 39.00
<b>17014453</b>	Multichannel All – 8 channels	108.40
<b>17014456</b>	Multichannel All – 12 channels	126.00
<b>30415513</b>	Liquidator 96 – 96 channels, includes box & round-trip UPS Ground shipping	1,850.00
<b>30456199</b>	BenchSmart 96 Head – 96 channels, includes round-trip UPS Ground shipping	1,850.00

ADVANCED PLAN – VOUCHER	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
When prior performance matters	✓	4 test weighings @ 10% and 100%	4 test weighings @ 10% and 100%	Standard Certificate

This plan includes all services provided under the Standard Plan.

Plus:

- As found calibration results – including failure confirmation

ORDER No.	DESCRIPTION	PRICE
<b>17014449</b>	Single channel	\$ 47.20
<b>17014454</b>	Multichannel All – 8 channels	146.60
<b>17014457</b>	Multichannel All – 12 channels	177.60
<b>30415514</b>	Liquidator 96 – 96 channels, includes box & round-trip UPS Ground shipping	2,200.00
<b>30456200</b>	BenchSmart 96 Head – 96 channels, includes round-trip UPS Ground shipping	2,200.00

ACCREDITED PLAN – VOUCHER	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
For complete ISO 17025 documentation	✓	4 test weighings @ 10%, 50% and 100%	10 test weighings @ 10%, 50% and 100%	ISO 17025 Certificate

This plan includes all services provided under the Advanced Plan.

Plus:

- As found NVLAP calibration results and failure confirmation
- ISO 17025 calibration certificate incl. uncertainty of measurement (certified by NVLAP), historical trending and Calibration label

ORDER No.	DESCRIPTION	PRICE
<b>17014450</b>	Single channel	\$ 91.00
<b>17014455</b>	Multichannel All – 8 channels	273.80
<b>17014458</b>	Multichannel All – 12 channels	306.40
<b>30415517</b>	Liquidator 96 – 96 channels, includes box & round-trip UPS Ground shipping	2,900.00
<b>30456201</b>	BenchSmart 96 Head – 96 channels, includes round-trip UPS Ground shipping	2,900.00

OPTIONAL RUSH SERVICE	GUARANTEED 48-HOUR IN-LAB TURNAROUND TIME WITH ROUND-TRIP UPS OVERNIGHT SHIPPING		
ORDER No.	DESCRIPTION	PRICE	
<b>30415518</b>	Rush Service Fee: Liquidator 96	\$ 525.00	
<b>17800408</b>	Rush Service Fee: BenchSmart 96 Head	525.00	

For more information on service selection,

please visit [www.RaininService.com](http://www.RaininService.com)

Phone orders: **1 800 4 RAININ (1 800 472 4646)**

Online orders: [www.ShopRainin.com](http://www.ShopRainin.com)